



clinell®

www.clinell.com

THE MOST CLINICALLY PROVEN DISINFECTANT WIPE IN THE WORLD

Patented formula

Clinell Sporidical Wipes come in a stable form that delivers peracetic acid on demand. Peracetic acid is a safe alternative to chlorine and is proven to provide similar or better performance against spores than chlorine¹.

Simple and easy to use

Just wet with water to activate the dry wipes. No more dilution errors. They do not generate toxic fumes and can be safely used next to patients.

The most powerful wipe in the world

Clinell Sporidical Wipes kill most known pathogenic microorganisms with a non-selective action. It has a proven 6 log spore kill, in 1 minute under dirty conditions². Proven to outperform other leading Sporidical wipes³, based on kill rates and contact times.

Contains detergent

Clinell Sporidical Wipes clean and disinfect. The peracetic acid remains effective in highly soiled conditions unlike chlorine based products which are inhibited by dirt and organic matter.

Safe to use on most surfaces

Clinell Sporidical Wipes are CE certified class IIa medical devices. They can be used on most non-invasive medical devices and equipment without corrosion, unlike chlorine solutions.

Cost efficient

Proven to save costs of £660,000 per annum⁴.

Proven to be faster

Proven to be faster and more effective at reducing spore counts than chlorine solutions⁵.

WHY CHOOSE CLINELL SPORICIDAL?

Improved cleaning and disinfection of room surfaces decreases the risk of healthcare associated infections²⁸.



Effective sporicidal agents

Enhanced environmental cleaning with sporicidal agents of rooms housing *Clostridium difficile* infected patients is warranted¹⁷. Key measures to prevent *C. difficile* transmission include correct cleaning and disinfection of the surfaces in hospital rooms daily and at discharge using a sporicidal disinfectant. For effective disinfection of *C. difficile*, a sporicidal product plus correct practices are essential¹⁸.

QACs ARE NOT SPORICIDAL

The number of chemical agents that possess sporicidal activity is limited to alkylating agents and oxidizing agents, such as peracetic acid, hydrogen peroxide and chlorine, the latter groups display more rapid sporicidal activity¹⁹. Quaternary Ammonium Compounds (QACs) are not sporicidal and are termed sporistatic^{20,21,22,23} which means they only inhibit spore germination and/or outgrowth. QACs do not kill spores.

CHLORINE TABLETS ARE INEFFICIENT

Chlorine solutions, made from powders or tablets break down over time. It is common practice to leave the prepared solution for many hours before using it, therefore at the time of use the solution can be ineffective.

Dilution errors can easily occur when preparing chlorine solutions. An error of using too much water when mixing up a solution can result in a disinfectant that is too weak and ineffective.

Conversely, using too little water results in a solution that is too strong, toxic and harmful to materials as well as to the user.

TOXICOLOGICAL HAZARDS


All disinfecting procedures must include a risk assessment of potential toxicological hazards²⁴.

Chlorine solutions, made from tablets and powders, are very toxic to both the user and the patients, emitting toxic, carcinogenic fumes and by-products. Excessive chlorine use has been shown to cause obstructive lung disease²⁵, shortness of breath, eye irritation, nasal complaints, cough and skin complaints²⁶.


HEALTH RISKS



Toxic fumes from chlorine based solutions can cause the following:



Lung damage, including cancer.



Eye irritation, nose and throat damage-causing coughs and shortness of breath.

Chronic inhalation of chlorine can increase the **risk of lung cancer**²⁷

CLINELL SPORICIDAL WIPES

Proven to reduce *Clostridium difficile* associated disease by **72%**⁴

Peracetic acid generating wipes activated by water, for surface disinfection and cleaning of non-invasive medical devices.



SPORICIDAL WIPES

Clinell Sporidical Wipes are a high level disinfectant wipe used specifically to target *Clostridium difficile* spores. It cleans and disinfects, providing a direct replacement and safe alternative to chlorine products^{1,5}.

Containing patented technology, Clinell Sporidical Wipes are designed for use on all surfaces of non-invasive medical devices. They are inactive when dry and with the addition of water they generate peracetic acid levels that are proven to kill most known microorganisms.

Kills at least
99.9999%
of spores in one minute,
in dirty conditions²

ANTIMICROBIAL ACTIVITY


Powerful disinfecting composition using a pH optimised peracetic acid and hydrogen peroxide combination, generated from sodium percarbonate and tetra acetyl ethylene diamine.

High power oxidative kill against all microorganisms including non-enveloped viruses and bacterial endospores.


Greater than 6 log kill (>99.9999%) of spores in 1 minute in dirty conditions². Conforms to EN1275, EN1276, EN13704, EN14348, EN14476, EN14561, EN14562 and EN14563.

CE 0473 Class IIa Medical Device

BACTERIA	TEST
<i>Acinetobacter baumannii</i>	EN13727
<i>Escherichia coli</i>	EN13727
<i>Pseudomonas aeruginosa</i>	EN13727
	EN14561
<i>Staphylococcus aureus</i>	EN13727
	EN14561
<i>Enterococcus hirae</i>	EN13727
	EN14561
<i>Klebsiella pneumoniae</i> (ESBL)	EN13727
<i>Enterococcus faecalis</i>	EN13727
<i>Enterococcus faecium</i> (VRE)	EN13727
SPORES	
<i>Clostridium difficile</i>	ASTM
	E2362-09
	EN13704
	Bab et al
<i>Bacillus subtilis</i>	
MYCOBACTERIA	
<i>Mycobacterium terrae</i>	EN14563
	EN14348
FUNGI	
<i>Candida albicans</i>	EN14562
	EN13727
<i>Aspergillus niger</i>	EN14562
	EN1275
VIRUSES	
Adenovirus	EN14476
Poliovirus	EN14476
Canine Parvovirus	EN14476
MERS CoV	EN14476



Dry before activation, each wipe consists of 2 non-woven fabrics bonded together with a unique powder combination in the middle layer.



Water activates the powder to generate high levels of hydrogen peroxide and peracetic acid.



This unique method of application ensures consistent efficacy, with no dilution errors.

PRODUCT INFORMATION	UNIT OF ISSUE	ORDER CODE	NHSSC
Sporidical Wipes	Pack of 25	CS25	VJT113
Sporidical Wall Mounted Dispenser for CS25		CS25D	-

USE DISINFECTANTS SAFELY. ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE. ALWAYS FOLLOW MEDICAL EQUIPMENT MANUFACTURERS CLEANING PROCEDURES AND GUIDELINES.

Clinell's industry leading aftercare and training services

Including a revolutionary training kit, specialised nurse trainers, educational materials, customised posters and dispenser installation.

Contact us on: 020 7993 0030 or info@clinell.com.

REFERENCES

1. Humphreys P, A systematic evaluation of a peracetic-acid-based high performance disinfectant. *Journal of Infection Prevention*. February, 2013, doi: 10.1177/1757177413476125.
2. Sporidical efficacy test, Hospital Infection Research Laboratory. May 2007. http://clinell.com/sites/default/files/data/Clinell_Sporidical_-_Efficacy_Test.pdf.
3. Siani H, et al. Efficacy of "sporidical" wipes against *Clostridium difficile*. *Am J Infect Control*. 2011 Apr;39(3):212-8. doi: 10.1016/j.ajic.2011.01.006.
4. Carter Y, Barry D. Tackling *C difficile* with environmental cleaning. *Nurs Times*. 2011 Sep 13-19;107(36):22-5.
5. Doan et al. Clinical and cost effectiveness of eight disinfection methods for terminal disinfection of hospital isolation rooms contaminated with *Clostridium difficile* 027. *J Hosp Infect*. 2012 Aug 14.
6. Sattar, S.A. and J.-Y. Maillard, The crucial role of wiping in decontamination of high-touch environmental surfaces: Review of current status and directions for the future. *Am J Infect Control*, 2013. 41(5, Supplement): p. S97-S104.
7. Panousi M et al. Evaluation of alcohol wipes used during aseptic manufacturing. *Lett Appl Microbiol*. 2009 May;48(5):648-51. doi: 10.1111/j.1472-765X.2009.02574.x. Epub 2009 Feb 18.
8. Maillard JY et al Emerging bacterial resistance following biocide exposure: should we be concerned? *Chemica Oggi*. 2009;27(3):26-8.
9. Karatzas KA, et al. Phenotypic and proteomic characterization of multiply antibiotic-resistant variants of *Salmonella enterica* serovar Typhimurium selected following exposure to disinfectants. *Appl Environ Microbiol*. 2008 Mar;74(5):1508-16. DOI: 10.1128/AEM.01931-07.
10. Gebel J et al. The role of surface disinfection in infection prevention. *GMS Hyg Infect Control*, Vol. 8(1), ISSN 1863-5245.
11. Meyer B, Cookson B. Does microbial resistance or adaptation to biocides create a hazard in infection prevention and control? *J Hosp Infect*. 2010 Nov;76(3):200-5. DOI: 10.1016/j.jhin.2010.05.020.
12. Maillard JY et al Emerging bacterial resistance following biocide exposure: should we be concerned? *Chemica Oggi*. 2009;27(3):26-8.
13. Bloss R et al. Adsorption of active ingredients of surface disinfectants depends on the type of fabric used for surface treatment. *J Hosp Infect*. 75 (2010) 56-61.
14. Engelbrecht, K., et al., Decreased activity of commercially available disinfectants containing quaternary ammonium compounds when exposed to cotton towels. *Am J Infect Control*, 2013, Epub ahead of print(0).
15. Goldsmith MT, Lattief MA, Friedl JL, Stuart LS. Adsorption of available chlorine and quaternary by cotton and wool Fabrics from disinfecting solutions. *Appl Microbiol* 1954;2:360-4.
16. Wren MW, Rollins MS, Jeanes A, Hall TJ, Coe'n PG, Gant VA. Removing bacteria from hospital surfaces: a laboratory comparison of ultramicrofibre and standard cloths. *J Hosp Infect* 2008;70:265e71.
17. Weber et al. Role of the environment in the transmission of *Clostridium difficile* in health care facilities. *Am J Infect Control*. 2013 May;41(5 Suppl):S105-10. doi: 10.1016/j.ajic.2012.12.009.
18. Cadnum et al. Transfer of *Clostridium difficile* Spores by Nonsporidical Wipes and Improperly Used Hypochlorite Wipes: Practice + Product = Perfection. *Infect Control Hosp Epidemiol*. 2013 Apr;34(4):441-2. doi: 10.1086/669871.
19. Maillard, J.Y., Innate resistance to sporicides and potential failure to decontaminate. *Journal of Hospital Infection*, 2011. 77(3): p. 204-209.
20. McDonnell, G., Antisepsis, disinfection, sterilization - types, action and resistance. 2007, Washington, D.C.: ASM Press.
21. Russell, A.D., Assessment of sporicidal efficacy. *International Biodeterioration & Biodegradation*, 1998. 41(3-4): p. 281-287.
22. Russell, A., Bacterial spores and chemical sporicidal agents. *Clinical Microbiology Reviews*, 1990. 3(2): p.99-119.
23. Siani, H., C. Cooper, and J.-Y. Maillard, Efficacy of "sporidical" wipes against *Clostridium difficile*. *American Journal of Infection Control*, 2011. 39(3): p. 212-218.
24. Gebel J et al. The role of surface disinfection in infection prevention. *GMS Hyg Infect Control*. 2013, Vol. 8(1), ISSN 1863-5245.
25. Sastre J et al. Airway response to chlorine inhalation (bleach) among cleaning workers with and without bronchial hyperresponsiveness. *Am J Ind Med*. 2011 Apr;54(4):293-9. doi: 10.1002/ajim.20912. Epub 2010 Oct 18.
26. LoVecchio F et al. Outcomes of chlorine exposure: a 5-year poison center experience in 598 patients. *Eur J Emerg Med*: June 2005 - Volume 12 - Issue 3 - pp 109-110.
27. Jappinen P et al. Cancer incidence among pulp and paper workers exposed to organic chlorinated compounds formed during chlorine pulp bleaching. *Scand J Work Environ Health*. 1991 Oct;17(5):356-9.
28. Weber et al. The role of the surface environment in healthcare-associated infections. *Curr Opin Infect Dis*. 2013 Aug;26(4):338-44. doi: 10.1097/QCO.0b013e3283630f04.

